Disposable Filter Needle

Approved by China NMPA, U.S. FDA and CE MDR

Version 3.1_2025-11-20



- *A built-in filter membrane with filtration efficiency up to 95% for insoluble particulates **
- Multiple needle lengths available to fit different vial sizes.
- Blunt needle type available for complete aspiration of medication liquid.
- Easy-to-use: no need to invert the vial when aspirating the drug.
- Hard PP shell and cap for primary packaging suitable for automatic assembling line.
- EB sterilization, no concern for ethylene oxide (EO) residue.

Procedures



Connect the filter needle to the Luer lock of the syringe.



Do not invert the vial, insert the filter needle into the vial and directly aspirate the liquid.



After aspirating the drug, cover the needle with the hard shell, and discard the needle.

Model & Specifications

Models	Gauge	Cannula Length	Pore Size of Filter	Tip Type	Primary Packaging
FN-B1838RW	18G	38mm	5μm	Blunt Regular Wall	PP Shell and Cap
FN-B1838RW15			15μm		
FN-B1838RW20			20μm		
FN-B1838RW25			25μm		
FN-B1850RW	18G	50mm	5μm		
FN-B1850RW15			15μm		
FN-B1850RW20			20μm		
FN-B1850RW25			25μm		
FN-L1938TW	19G	38mm	5μm	Lancet Thin Wall	
FN-L1938TW15			15μm		
FN-L1938TW20			20μm		
FN-L1938TW25			25μm		
FN-L1950TW	19G	50mm	5μm		
FN-L1950TW15			15μm		
FN-L1950TW20			20μm		
FN-L1950TW25			25μm		

* Data from our laboratory.

NEXTECH™

Hangzhou Qiantang-Longyue Biotechnology Co., Ltd.
Website: www.nextech-x.com/en
Email: sales1@nextech-x.com
Tel: 0086-0571-89150121, 0086-18017582758