Vial Adapter


Vented Type & Drug-Filtered Type


FDA 510(k) clearance and CE MDR certification



Vented Type

- One piercing spike with two lumens:
  - One for transfer of liquid medication between the drug vial and a syringe
  - Second for enabling pressure equilibrium of the vial by introducing air that has passed through a 0.2μm hydrophobic air filter into the vial
- Options for siliconized spike
- Available for vial-first or syringe-first packaging





Drug-Filtered Type

- One drug filter built in the liquid lumen of the Vial Adapter Basic Type
- Available in 5μm, 15μm and 25μm pore sizes
- Effectively removes insoluble particulates for injection safety

Version 1_2026-08-04

NEXTECH™

Hangzhou Qiantang Longyue Biotechnology Co., Ltd.
Rooms 302 and 104, Building 12, Block 1, No.619 Wangmei Road,
Linping Street, Linping District, Hangzhou City
Website: www.nextech-x.com/en
Tel: 0086-0571-89150121, 0086-18017582758