Vial Adapter *Multi-Dose*

Approved by Australian ARTG, Health Canada and Saudi FDA.

Version 1_2026-07-20



Features and Uses

- Consists of the vial adapter body with a single piercing spike and a swabable valve.
- Intended to allow multiple needle-less accesses to an injection medication vial for the purpose of facilitating the withdrawal or addition of drugs/solutions from or to a vial.



NEXTECH™

I Product Series

❶ Basic Vial Adapter *Multi-Dose*

- One piercing spike with a lumen for transfer of drug/solution between the vial and a syringe.



❷ Vented Vial Adapter *Multi-Dose*

One piercing spike with two lumens:

- One for transfer of drug between the vial and a syringe
- Second for automatic pressure equilibrium of the vial by introducing air that passes through a 0.2μm hydrophobic air filter into the vial.



❸ Drug-Filtered Vial Adapter *Multi-Dose*

- One piercing spike with a lumen for transfer of drug between the vial and a syringe
- One in-line drug filter to eliminate insoluble particles from the solution



Certificate Information

- Australian ARTG Certificate DV-2025-MC-14825-1
- Health Canada Medical Device Licence 113412
- Saudi FDA Medical Device Certificate ME0000000773SFDA00010

NEXTECH™

Hangzhou Qiantang Longyue Biotechnology Co., Ltd.
Rooms 302 and 104, Building 12, Block 1, No.619 Wangmei Road,
Linping Street, Linping District, Hangzhou City
Website: www.nextech-x.com/en
Tel: 0086-0571-89150121, 0086-18017582758